



Tabletop Refraction System TS-310







Triple-performance combination

Operation-oriented control box

Intuitive color LCD touch screen

The 5.7-inch color LCD touch screen displays all data with high visibility. Simple, understandable, and comfortable user interfaces ensure an effective operation.



Dial button with S / C / A mode key

Ergonomically focused, the main control dial, which has S / C / A mode key in its center, allows for quick and smooth refractive value changes.



Built-in high-speed line printer

The control box includes a built-in, high-speed printer which automatically outputs all measurement data in a format that provides easy interpretation and explanation to patients.



Sophisticated refractor

Ergonomic design

Beautiful ergonomic design enhances a stress-free examination environment for both operator and patient, while maintaining superior accuracy.



Smooth lens change

Speedy, smooth and quiet lens changing enables operator to make minute changes of prescription at the touch of a button.



Antifogging protective glass, detachable forehead rest and face shields, support easy maintenance.



Reliable chart unit

Same charts for far and near vision measurement

The TS-310 uniquely uses the same high resolution charts for both far and near testing. Switching between far and near is smoothly achieved with a single button push.

High precision, large, color LCD

The brilliant clarity of the LCD, allows visual acuity at 5 m and near visual acuity at 40 cm to be measured with the same accuracy as actual distances.

Contrast test


To support patients who have undergone cataract or refractive surgery, the TS-310 has the ability to measure contrast sensitivity at three different levels below the normal threshold.



Black and white inversion function

The VA chart can be selected from black on white to white on black.





TS-310 Chart Type*

	Туре Т	
	Distance / near	Near
Letter	0.04 to 1.5	0.1 to 1.0
Number	0.04 to 1.5	0.1 to 1.0
Tumbling E	0.1 to 1.5	0.1 to 1.0
Landolt C	_	_
Children	0.1 to 1.0	_
ETDRS style letters	0.32 to 2.0	_
Letters and Numbers	_	0.63
Astigmatism clock dial		_
Dots	-	_
Red-green	8 🖷 8	_
Cross grid		_
Binocular balance	5 > * * *	_
Duochrome balance	60 09	_
Phoria		_
	+ +	
Phoria with fixation point		_
Vertical line	_	
Horizontal line	_	
Vertical coincidence		_
Horizontal coincidence	_	_
Schober	((+))	_
Stereo	-11	_
Worth four dot	+,+	_
Fixation point	•	

^{*} In addition to the Type T chart, other chart types are available.

User-friendly features

Symmetry design

Clean symmetry design, without a refractor pole, allows examinations from both right and left side of patients. In this way, multiple TS-310 units can be positioned with the greatest flexibility.

Flexible measurement style

Regardless of personal measurement style, or sitting / standing ergonomic preferences, the TS-310 ensures comfortable examinations.






Easy installation

The TS-310 does not require complex angle and distance adjustments during installation, which can be achieved by any staff member without expert knowledge. It can be easily transported and installed in any facility including optometric offices, mobile optometric vehicles, and care facilities for temporary use.

Operational efficiency

The TS-310 embodies numerous ingenuities, which support both accuracy and functional flexibility, to create a new standard in examination efficiency.



- 1. Pilot lamp for power on / off / auto off mode
- 2. Power switch
- 3. Control box tray (optional)


- 4. VD check adapter (standard / additional adapters are optional.)
- 5. Refractor head vertical motion switch
- 6. Pupil position check window / cover

Scalable according to your needs

Due to the compact and simple body shape, the TS-310 fits any table or location.


NIDEK refraction products allow for quick and easy wireless data transfer* using the Eye Care card, Bluetooth, WLAN or infrared communication. This is helpful for making a simple refraction system without complicated wired connection.

*The specifications for wireless data transfer differ according to each product and from country to country. The requirements also differ depending on the method of wireless data transfer.



Compact footprint

Super space-saving design of the TS-310 (from 0.4 to 0.7 m wide, 0.5 m deep), offers flexible room arrangement without compromising the precision of conventional refractions. The compact footprint is unchanged, whether the refractor head is in the raised or lowered position.



TS-310 Specifications

Chart type T, UK, M, P, PhM Refraction distance Distance: 5 m, Near: 40 cm Mask Horizontal line, Vertical line, Single letter Filter Red (Green) Binocular vision test Red-green, Prism Refractor Measurable range Sphere -19.00 to +16.75 D (0.25 / 0.5 to 3 D increments) Cylinder 0.00 to ±6.00 D (0.25 / 1 D increments) Axis 0 to 180° (1/5° increments) PD 48 to 80 mm Fish on Mark (Green) 40 to 80 mm (100% convergence possible PD) Prism 0 to 20∆ (0.1 / 0.5 / 2∆ increments) Auxiliary lens 40.25 D Cross cylinder lens 40.25 D Octuder Available Pinhole plate Red (Green filter Red ross cylinder lens Available Sherical lenses for retinoscope Red (Ersen filter) Red maddox rod Right eye: horizontal. Left eye: yerical Dissociation prism Right eye: 6a8U / Left eye: 10.38, Right eye: 3.3BD / Left eye: 3.3BU Visual field 40° VD = 12 mm), 39° (VD = 13.75 mm) Vertex distance marking 2.1, 13.75, 16,	- S To Specifications	
Refraction distance Distance: 5 m, Near: 40 cm	Chart	
Horizontal line, Vertical line, Single letter Red / Green Red / G	Chart type	T, UK, M, P, PhM
Filter Binocular vision test Red / Green Red-green, Prism Refractor	Refraction distance	Distance: 5 m, Near: 40 cm
Binocular vision test Red-green, Prism Refractor Measurable range Sphere -19.00 to +16.75 D (0.25 / 1.0 increments) Cylinder 0.00 to ±6.00 D (0.25 / 1.0 increments) Axis 0 to 180° (17/5° increments) PD 48 to 80 mm Auxiliary lens ±0.25 D Cross oflinder lens ±0.25 D Occluder Available Pinhole plate a0 zmm Red / Green filter Right eye: red, Left eye: green PD check lens ±0.50 D Spherical lenses for retinoscope 45.7 + 2.0 D Red maddox rod Right eye: 68.80 / Left eye: vertical Q bissociation prism Right eye: 68.80 / Left eye: 10.8B, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment ±2.5° Refractor arm Lectrically-driven Up-and-down of refractor arm 190 mm Up-and-down of refractor arm 190 mm Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x × 87; 10) (W) x 20 (W) x 70 (H) mm 17.6 (W) x 20 (W) x 3	Mask	Horizontal line, Vertical line, Single letter
Refractor Measurable range Sphere -19.00 to ±16.75 D (0.25 / 0.5 to 3 D increments)	Filter	Red / Green
Measurable range Sphere	Binocular vision test	Red-green, Prism
Sphere	Refractor	
Cylinder	Measurable range	
Cylinder	Sphere	-19.00 to +16.75 D (0.25 / 0.5 to 3 D increments)
Axis	Cylinder	0.00 to ±6.00 D (0.25 / 1 D increments)
PD		0 to 180° (1°/5° increments)
Frism Auxillary lens Cross cylinder lens Occluder Available Pinhole plate Po check lens Fixed cross cylinder lens Spherical lenses for retinoscope Red maddox rod Dissociation prism Visual field Vertex distance marking Level adjustment Vertex distance marking Level adjustment Refractor arm Up-and-down of refractor arm Up-and-down of refractor arm Power supply Ac 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 Auxilable ±0.50 D +1.5 / +2.0 D		
Prism		
Auxiliary lens Cross cylinder lens Occluder Available Pinhole plate Red / Green filter PD check lens Fixed cross cylinder lens Spherical lenses for retinoscope Red maddox rod Dissociation prism Right eye: horizontal, Left eye: vertical Dissociation prism Right eye: 6aBU / Left eye: 10aBl, Right eye: 3ABD / Left eye: 3ABU Visual field Vivual field Forehead rest adjustment Vertex distance marking Level adjustment Power supply Spherical printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)* Mass*2 Mass*2 Mass*2 Mass*2 Masi*But eye: cab Left eye: green Available a valiable ver red, Left eye: green Available ver ted ver sab Die left eye: oabple of the set of t	Prism	
Cross cylinder lens Occluder Occluder Available Pinhole plate Red / Green filter PD check lens Fixed cross cylinder lens Spherical lenses for retinoscope Red maddox rod Dissociation prism Visual field Forehead rest adjustment Level adjustment Level adjustment Level adjustment Level adjustment Level adjustment Up-and-down of refractor arm Up-and-down of refractor arm Display Printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 20 x (D) x 35.7 (H)" Mass*2 Mass*2 Mass*2 Mailable Available Ava		V to 204 (c. 17 0.07 24 interments)
Occluder Pinhole plate Pinhole plate Red / Green filter Red / Green filter PD check lens Fixed cross cylinder lens Spherical lenses for retinoscope Red maddox rod Dissociation prism Right eye: horizontal, Left eye: vertical Poissociation prism Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment Vertex distance marking 12, 13.75, 16, 18, 20 mm Level adjustment Level adjustment Up-and-down of refractor arm In high speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.	•	+0.25 D
Pinhole plate Red / Green filter PD check lens PD check lens Fixed cross cylinder lens Spherical lenses for retinoscope Red maddox rod Dissociation prism Visual field Vorument Vertex distance marking Level adjustment Level adjustment Vertex distance marking Level adjustment Up-and-down of refractor arm Up-and-down of refractor arm Up-and-down of refractor arm Pinter Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 Mass*2 Mass*2 O 2 mm Right eye: red, Left eye: green Available 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment 4 ± 2 mm 12, 13.75, 16, 18, 20 mm 4 ± 2.5° Refractor arm Electrically-driven 190 mm Control box Display Printer High speed line printer Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" 32 kg 70.5 lbs.		
Red / Green filter		
PD check lens Fixed cross cylinder lens 5pherical lenses for retinoscope Red maddox rod Dissociation prism Visual field Forehead rest adjustment Vertex distance marking Level adjustment Up-and-down of refractor arm Up-and-down of refractor arm Display Printer Power supply Power consumption Dimensions*1 Mass*2 Mass*2 Available 4.0.5 0 D Available 4.0.5 0 D Available 4.0.5 0 D Right eye: 40.5 0 D Right eye: vertical PoliApl, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBU / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.5 0 D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.0 N D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.0 N D Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.0 N D Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 4.0.0 N D Right eye: 6ΔBU / Left eye: 3ΔBU	•	
Fixed cross cylinder lens Spherical lenses for retinoscope Red maddox rod Dissociation prism Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment 14 ±2 mm Vertex distance marking Level adjustment 2.15°, 16, 18, 20 mm Level adjustment Up-and-down of refractor arm Up-and-down of refractor arm Up-and-down of refractor arm Display Printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)* The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)* 32 kg 70.5 lbs.		
Spherical lenses for retinoscope Red maddox rod Dissociation prism Right eye: 6ABU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment Vertex distance marking 12, 13.75, 16, 18, 20 mm Level adjustment 42.5° Electrically-driven Up-and-down of refractor arm Up-and-down of refractor arm 190 mm Control box Display 5.7-inch color Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H) mm 17.6 (W) x		
Red maddox rod Dissociation prism Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment 14 ±2 mm Vertex distance marking 12, 13.75, 16, 18, 20 mm Level adjustment 42.5° Refractor arm Level adjustment 42.5° Refractor arm Up-and-down of refractor arm Up-and-down of refractor arm Display Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" 32 kg 70.5 lbs.		
Dissociation prism Right eye: 6ΔBU / Left eye: 10ΔBI, Right eye: 3ΔBD / Left eye: 3ΔBU Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment 14 ±2 mm Vertex distance marking 12, 13.75, 16, 18, 20 mm Level adjustment ±2.5° Refractor arm Electrically-driven Up-and-down of refractor arm 190 mm Control box 5.7-inch color Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)* The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" 32 kg 70.5 lbs. 70.5 lbs.	·	
Visual field 40° (VD = 12 mm), 39° (VD = 13.75 mm) Forehead rest adjustment 14 ±2 mm Vertex distance marking 12, 13.75, 16, 18, 20 mm Level adjustment ±2.5° Refractor arm Electrically-driven Up-and-down of refractor arm 190 mm Control box 5.7-inch color Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.		
Forehead rest adjustment Vertex distance marking Level adjustment Level adjustment Refractor arm Up-and-down of refractor arm Up-and-down of refractor arm Display Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 Mass*2 Mass*2 Mass*2	·	
Vertex distance marking 12, 13.75, 16, 18, 20 mm Level adjustment ±2.5° Refractor arm Electrically-driven Up-and-down of refractor arm 190 mm Control box Display Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.		
Level adjustment Refractor arm Electrically-driven 190 mm Control box Display Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 30.0 (H)" The refractor head vertion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 Mass*2 Mass*2	•	14 ±2 mm
Refractor arm Up-and-down of refractor arm 190 mm Control box Display Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 20.4 (D) x 35.7 (H) " Mass*2 Mass*2 32 kg 70.5 lbs.	Vertex distance marking	12, 13.75, 16, 18, 20 mm
Up-and-down of refractor arm 190 mm Control box Display Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 Mass*2 32 kg 70.5 lbs.	Level adjustment	±2.5°
Control box Display Printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.	Refractor arm	Electrically-driven
Display Printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 Mass*2	Up-and-down of refractor arm	190 mm
Printer High speed line printer Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom.	Control box	
Power supply AC 100 to 240 V, 50 / 60 Hz Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H) " The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H) " Mass*2 32 kg 70.5 lbs.	Display	5.7-inch color
Power consumption 130 VA Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H) " The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H) " Mass*2 32 kg 70.5 lbs.	Printer	High speed line printer
Dimensions*1 The refractor head vertical motion unit is at the bottom. 446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H) " The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H) " Mass*2 32 kg 70.5 lbs.	Power supply	AC 100 to 240 V, 50 / 60 Hz
446 (W) x 519 (D) x 763 (H) mm 17.6 (W) x 20.4 (D) x 30.0 (H) " The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H) " Mass*2 32 kg 70.5 lbs.	Power consumption	130 VA
17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.	Dimensions*1	The refractor head vertical motion unit is at the bottom.
17.6 (W) x 20.4 (D) x 30.0 (H)" The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.		446 (W) x 519 (D) x 763 (H) mm
The refractor head vertical motion unit is at the top. 446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.		
446 (W) x 487 (D) x 907 (H) mm 17.6 (W) x 19.2 (D) x 35.7 (H)" Mass*2 32 kg 70.5 lbs.		
17.6 (W) x 19.2 (D) x 35.7 (H)" 32 kg 70.5 lbs.		·
Mass*2 32 kg 70.5 lbs.		
70.5 lbs.	Mass*2	
	101033 -	g .
	Standard accessories	Forehead rest, Face shields, VD check adapter, Stylus pen, Printer paper, Dust cover, Power cord, Cover plate, Cap, Screw,
Hexagonal wrench, Phillips screwdriver	Standard accessories	
	Ontional accessories	Eye Care card, Communication cable for an auto refractometer, Communication cable for a computer, Control box tray
Lye Care card, Communication cable for an auto retractometer, Communication Cable for a Computer, Control box to	Optional accessories	Lye care card, Communication cable for an auto refractioneter, Communication cable for a Computer, Control box tray

^{*1} Control box excluded.

Product / Model name: Tabletop Refraction System TS-310

All LCD images are simulated.

Brochure and listed features of the device are intended for non-US practitioners.

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.





800.874.5274 www.marco.com 11825 Central Parkway Jacksonville, FL 32224 USA

HEAD OFFICE (International Div.) 34-14 Maehama, Hiroishi Gamagori, Aichi 443-0038, JAPAN

TEL: +81-533-67-8895 URL: http://www.nidek.com [Manufacturer]

TOKYO OFFICE (International Div.)

3F Sumitomo Fudosan Hongo Bldg., 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN TEL: +81-3-5844-2641 URL: http://www.nidek.com

NIDEK INC. 47651 Westinghouse Drive, Fremont, CA 94539, U.S.A.

TEL: +1-510-226-5700 +1-800-223-9044 (US only) URL: http://usa.nidek.com NIDEK S.A.

Europarc, 13 rue Auguste Perret, 94042 Créteil, FRANCE TEL: +33-1-49 80 97 97 URL: http://www.nidek.fr NIDEK TECHNOLOGIES S.R.L. NIDEK (SHANGHAI) CO., LTD. Via dell'Artigianato, 6/A, 35020 Albignasego (Padova),

TEL: +39 049 8629200/8626399 URL: http://www.nidektechnologies.it

Rm3205,Shanghai Multi Media Park, No.1027 Chang Ning Rd, Chang Ning District, Shanghai, CHINA 200050, TEL: +86 021-5212-7942 URL: http://www.nidek-china.cn NIDEK SINGAPORE PTE. LTD. 51 Changi Business Park Central 2, #06-14, The Signature 486066, SINGAPORE TEL:+65 6588 0389

^{*2} Refractor head and control box included.