

KM-500 HAND-HELD AUTOMATIC KERATOMETER

*Objective, infrared
measurement means more
accurate, repeatable results.*



*Hand-held operation is
completely portable, making it
easier to measure children and
physically-challenged patients.*



*Magnifying window allows
the operator to stand at a
comfortable measuring
distance from the patient.*



*Extremely fast and simple,
the KM-500 automatically
takes the measurement
after reaching the proper
alignment distance.*



*At the touch of a button,
both millimeter and diopter
measurements are displayed,
as well as average K-readings
and corneal cylinder.*



Automatic keratometry is a snap with the **KM-500 HAND-HELD KERATOMETER**. The KM-500's unique viewing and alignment window lets you observe the targeted eye from a non-threatening distance. There is no need to dance "cheek-to-cheek" with your patient.

Measurement is simple. View the measuring eye through the window and move the KM-500 towards the patient. When the mire ring appears on the cornea, slowly move the unit closer. The mire ring will "flash" and the KM-500 will fire automatically when the proper distance is reached. Readings appear instantly on the LCD display. Both left and right readings display simultaneously.

The optional printer will allow you to print results by simply holding the KM-500 above the printer and transferring the data by infrared transmission. There is no need to dock the unit to the printer. Two rechargeable batteries provide continuous uninterrupted use through any busy schedule. Each battery provides 60 minutes of constant readings.



SPECIFICATIONS	
ALIGNMENT	Reticle ring
MEASUREMENT	Automatic
MEASURABLE RANGE	33–67.5 diopters
CYLINDER DISPLAY	Minus Cylinder Only
RADIUS DISPLAY	R1 Minor Meridian R2 Major Meridian
POWER TYPE	Battery, replaceable packs
POWER LIFE	Good for 1 hour continuous use
RECHARGE TIME	90 minutes
MODELS	Hand-held
WEIGHT	1.5 lbs.
OPTIONAL PRINTER	Wireless
DISPLAY	LCD, R & L eye

PLAY IT SAFE.
Protect your warranty. Buy your Marco product only from an authorized Marco distributor. Call Marco for the name of the distributor nearest you.



Marco-Certified Quality assures you of the very finest vision products, offering advanced technology, industry-leading durability and enhanced productivity for your practice. Marco-Certified Quality is your mark of confidence.

SC05000 BR-KM500

T H E L E A D E R I N V I S I O N D I A G N O S T I C S ®



11825 Central Parkway | Jacksonville, FL 32224 | 800 874-5274 | 904 642-9330 | www.marco.com