





AT-9

**APPLANATION
TONOMETER**

NEW



MARCO AT-9 TONOMETER SPECIFICATIONS

Specifications	
Slit Lamp Model	Ultra M3, Ultra M5
Tonometer Mounts	M3 #1688, M5 #1687
Measurement method	Applanation type (observed with left eye on the slit lamp)
Depression Area	7.354mm ²
Diameter of Measurement Depression Area	3.06mm
Measurement Range	0 ~ 80mmHg
Minimal Scale	2mmHg
Equation for Conversion from mmHg to kpa	1mmHg = 101325/760Pa = 133.3Pa 100Pa = 0.1kPa


Check Weight & Check Weight Holder



Two Measuring Prisms Standard



AT-9 Tonometer



Smoother rotation and operation allows for more accurate positioning and approach at the patient cornea.

