

**Marco
SC-1600**

LCD ACUITY CHART SYSTEM



- High resolution LCD screen provides a wide range of popular Snellen, ETDRS, and stereo testing options for binocular vision examinations.
- A dedicated, solid-state refracting device with an extreme emphasis on precise visual acuities.
- White background with uniform brightness for greater contrast and sharper letter edge quality.
- Incorporates all of the advanced features found on the first generation SC-2000 LCD Acuity Chart System and the popular series of Marco automated chart projectors.
- Contrast sensitivity testing measures contrast at three distinctive levels below the normal threshold.
- Randomization feature helps to prevent patient letter/character memorization.
- User-friendly remote control includes an LCD chart display window for convenient operator reference.
- Ability to add up to 18 personalized images on removable SD card.




SC-1600 SPECIFICATIONS



LCD	
Type	17-inch SXGA LCD with long-lasting backlight
Resolution	1280 x 1024 pixels
Chart Window Size	13.6" (W) x 10.86" (H)
Brightness	200 cd/m2 (5 level variable)
Charts	Type M, Type UK, Type P, Type T
Distance of Installation	2.5-7m (8'-23') (in the Mirror mode; 5 to 7mm), 0.5m increments
Chart Selection Devices	Remote control or control panel for TRS (RT-5100) refraction system
Memory Card	CF card / Built-in card slot
Power Supply	AC100 V- 240 V, 50 / 60 Hz
Power Consumption	100 VA
Dimensions & Weight	20.5" (W) x 18" (H) x 1.9" (D) / 11 lbs.
Standard Accessories	Anchor Screws (x4), Mounting Plate, Tilting Mounting Bracket

SC-1600 Chart Types



All automated Marco technologies share data and integrate with leading EMR programs via Marco Connect software.

Slide Show, Contrast Sensitivity, ETDRS Charts

Pediatric Snellen Charts

Number Charts

Special Tests

Program Random Line Functions

ETDRS Detail Sticker


Chart Display Window

Illiterate E Snellen Charts

Alphabet Snellen Charts

Black/White Background Switch

Maneuverable Masking Options (Vertical, Horizontal, Single Character)



Marco-Certified Quality assures you of the very finest vision products, offering advanced technology, industry-leading durability and enhanced productivity for your practice. Marco-Certified Quality is your mark of confidence.



800.874.5274
www.marco.com

