ME-CHECK MEIBOMIAN GLAND IMAGING



Capture high-quality meibomian gland images in your exam room with ME-CHECK.

FEATURES:



HIGH-RESOLUTION
MEIBOMIAN GLAND IMAGING



EFFORTLESSLY MOUNTS TO MOST SLIT LAMPS



EASY, INTUITIVE IMAGING SOFTWARE



PLUGS DIRECTLY INTO YOUR PC USB PORT FOR SIMPLE EMR UPLOAD

Simply mount it on your slit lamp, align, click, and easily save to your EMR or imaging database.

ME-CHECK INCLUDES:

- High-resolution infrared camera with infrared illumination
- Universal pin mount for slit lamp attachment
- Storage dock
- · PC compatible imaging software
- 9' USB connection cable
- · Lid eversion tool



1.800.874.5274 | MARCO.COM