

# MARCO BRAVO CHAIR






# MARCO Bravo Chair

For years, practitioners have asked for a simple-featured/lower-cost exam chair for everyday primary patient care. The new BRAVO chair from Marco embraces this basic concept in a contemporary design that compliments our established family of classical products. By eliminating the recline mechanism and protruding “elbow” headrest, the BRAVO chair can be comfortably positioned against the wall, maximizing the refraction distance in small exam rooms. Convenient, pop-on cushions feature a high curved cushion back that provides a comfortable patient headrest position for eye drops or indirect exams. And with dozens of different exam lane configurations to choose from, you can easily customize your own economical Marco lane package.

- Simple-featured/lower cost exam chair designed for everyday primary patient care.
- Up/down operations are controlled by corded foot pedal or from numerous locations on instrument stand.
- High, curved cushion back provides a comfortable patient headrest position for eye drops or indirect exams.
- Modern, high-quality electronic lift mechanism ensures reliable, long-term elevation control.
- Increased elevation range easily accommodates all patients and practitioners.
- No recline function allows user to maximize refraction distance in small exam rooms.
- Solid, one-piece steel frame accommodates easy pop-on cushions.
- Scuff-resistant, high-impact covers on arm rests and chair base.
- Standard Marco ten-step, hand-finished painting process.



**Complete lane packages and sophisticated integration with other Marco products are also attractive options.**



**High, curved cushion back provides a comfortable patient headrest position. Scuff-resistant arm rests are also a standard feature.**



**Modern, high quality lift mechanism provides extended elevation range.**

Specifications	
Weight	185 lb
Stroke	15.75"
Lifting Capacity	Approx. 350 lbs.
Electrical	120v, 50/50 Hz