


VA-1 ACUITY SYSTEM

VA-1 LCD VISUAL ACUITY SYSTEM



COMPREHENSIVE VISUAL ACUITY SOLUTION

- Multiple optotype selections with letters, E's, numbers, Landolt rings, HOTV, Snellen, and children's VA charts
- All acuity slides presented with ETRS spacing
- Contrast sensitivity testing
- Crowding bars for pediatric testing
- R/G binocular tests
- Multimedia system: animations, cartoons, and education slides included and more can be added
- Specialty slides including color tests


ADDED FEATURES:

- Intuitive remote controller
- Randomization of optotypes
- Custom sequences can be programmed
- Memory button: allows instant slide selection of your selected slide
- High resolution 24" LCD screen: provides optimal brightness for clear presentation of slides

OPTOTYPES



BINOCULAR TESTS



EXTENDED ETDRS CHARTS



ASTIGMATISM CHARTS



CONTRAST SENSITIVITY





REVERSED



DUOCHROME TEST




24 COLOR TESTS



CROWDING BARS



MULTIMEDIA



VA-1 ACUITY SYSTEM SPECIFICATIONS

Display	24in. WUXGA TFT LCD (1920 x 1200px)
Brightness	Maximum 300cd/m2
Dimension & Weight	(W)568 x (H)372 x (D)54mm, Approximately 4.0kg
Power Supply	DC 12V 5A
Measuring Distance	9ft. to 26ft. (2.7 to 8 meters)
External Interface	4x USB, RS-232, 1x Ethernet (10/100) , 1x Serial

CONTACT MARCO TO LEARN MORE.

Marco Ophthalmic | (800) 874-5274