



BRAVO 2

THE NEXT GENERATION INSTRUMENT STAND

The Bravo 2 Instrument Stand from Marco combines a fresh, new updated design with our high standards of performance, function, and durability. With a strong emphasis on simplicity and economy, the Bravo 2 stand provides the perfect blend of everyday efficiency and value.

FEATURES:

- » Attractive, modern design compliments any practice/office environment
- » Simplified membrane console panel provides convenient exam chair and overhead lamp controls
- » Robust lower slit lamp and refractor arms provide smooth, secure movement of mounted instruments



The Leader in Advanced Technology



INSTRUMENT CONSOLE PANEL

A simplified, redesigned membrane console panel provides controls for exam chair elevation and auto return, overhead lamp with adjustable rheostat, and three rechargeable instrument wells that accommodate all brands of handheld instruments.



OVERHEAD LAMP

Consistent with all Marco stands, the overhead lamp includes a long-lasting LED bulb that provides bright, direct illumination for general lighting or focused minor procedures. Durable construction allows for stable, one-handed adjustments in any direction.



COUNTERBALANCED REFRACTOR ARM

Carefully built with a natural, comfortable motion, the counterbalanced refractor arm includes a solid locking/release lever that ensures precise adjustments and secures locking of the refractor head in all three dimensions.



COUNTERBALANCED SLIT LAMP ARM

Creatively designed and engineered, the robust counterbalanced slit lamp arm provides an expanded range of motion to better accommodate a wide variety of patients, including handicap patients. A conveniently located locking/release mechanism provides smooth up/down travel and locking of the slit lamp. An electronic outlet and USB charging port are also included.

BRAVO 2 INSTRUMENT STAND TECHNICAL SPECIFICATIONS

STAND

- » Base: 22" x 23.25"
- » Height Console: 33.5"
- » Height Pole: 68"

SLIT LAMP ARM

- » Capacity: Up to 60 lbs
- » Range: 13.5"

REFRACTOR ARM

- » Capacity: Up to 20 lbs
- » Range: 12.5"

Contact Marco First. MARCO.COM | 800.874.5274

Marco Ophthalmic, 11825 Central Parkway, Jacksonville, Florida 32224

A DIVISION OF  ADVANCING EYECARE