

# **MARCO CHAIRS & STANDS**

The Difference is Marco

# **MARCO E-Z TILT**

The Marco E-Z Tilt Chair combines the most simplified manual reclining motion with a new modern, sophisticated appearance. The Marco E-Z Tilt also incorporates the only double hydraulic cylinder system in the industry, creating an extended range of elevation, from extra low to extra high, helping to accommodate a full range of patient and practitioner sizes. Convenient controls for chair reclining, elevation, rotation, and headrest adjustments are fingertip and foot operated.





# **MARCO ENCORE AUTOMATIC**

The Marco Encore Automatic Chair features an exclusive hydraulic system that controls both the elevation and recline motion. Easy-to-use backlit fingertip switches control all functions of the chair, as well as programmable recline positions, autorecline/return, and safety functions. A corded foot pedal and single-handed operation headrest are also standard. Optional surgical accessories are available for minor surgical procedures.



EASY, "DIAL-IN"
PROGRAMMING CONTROL
WITH SAFETY SWITCH TO
DISABLE CHAIR OPERATION

# **MARCO ENCORE MANUAL**

Considered our standard, everyday exam chair, the Marco Encore Manual Chair incorporates a simple, pneumatic-assisted recline mechanism for easy patient adjustments. Convenient fingertip switches on the side of the chair instantly control chair elevation and auto return. A corded foot pedal and single-handed operation headrest are also standard.





# **TECHNICAL SPECIFICATIONS**

#### **EZ TILT EXAMINATION CHAIR**

- » Base: 22.25" x 23.85"
- » Lift Capacity: Approx 400 lbs
- Solution Seat Height: 20.5" low 33" high
- » Rotation: 270° degrees
- » Recline Mechanism: Tilting recline
- » Input Voltage: 120VAC 60HZ
- » Fuses: 8A 250V
- » Avaliable Colors: Black, Charcoal, Navy, Silver, Silver & Black two tone, Silver & Navy two tone

## **BRAVO 2 EXAMINATION CHAIR**

- » Base: 22" x 23.25"
- » Lift Capacity: Approx. 350 lbs
- » Ground to Seat Height: 21.5" low 30" high
- » Rotation: 270° degrees
- » Recline Mechanism: Tilting recline
- » Input Voltage: 100-120VAC 60HZ
- » Fuses: 8A
- » Avaliable Colors: Black, Charcoal, Navy

## **ENCORE MANUAL EXAMINATION CHAIR**

- » Base: 22.25" x 23.85"
- » Lift Capacity: Approx. 375 lbs
- Solution Seat Height: 21" low 32" high
- » Rotation: 270° degrees
- » Recline Mechanism: Manual Lay Back
- » Input Voltage: 120VAC 60HZ
- » Fuses: 8A 250V
- » Avaliable Colors: Black, Charcoal, Navy, Silver, Silver & Black two tone, Silver & Navy two tone

#### **DELUXE INSTRUMENT STAND**

#### UNIT

- » Base: 21.25" x 21.25"» Console Height: 33.5"
- » Power: AC 120V (50/60 Hz)

#### SLIT LAMP ARM

- » Capacity: Up to 60 lbs
- » Range: 13"

#### REFRACTOR ARM

- » Capacity: Up to 12 lbs
- » Range: 8"

## **BRAVO 2 INSTRUMENT STAND**

#### UNIT

- » Base: 22" x 23.25"
- » Height Console:33.5"
- » Height Pole: 68"

#### SLIT LAMP ARM

- » Capacity: Up to 60 lbs
- » Range:13.5"

#### REFRACTOR ARM

- » Capacity: Up to 20 lbs
- » Range:12.5"

### **ENCORE AUTOMATIC EXAMINATION CHAIR**

- » Base: 22.25" x 23.85"
- » Lift Capacity: Approx. 375 lbs
- Solution of the Ground to Seat Height: 21" low 32" high
- » Rotation: 270° degrees
- » Recline Mechanism: Auto Lay Back
- » Input Voltage: 120VAC 60HZ
- » Fuses: 8A 250V
- » Avaliable Colors: Black, Charcoal, Navy, Silver, Silver & Black two tone, Silver & Navy two tone

Contact Marco First. MARCO.COM | 1.800.874.5274

Marco Ophthalmic, 11825 Central Parkway, Jacksonville, Florida 32224

A DIVISION OF ADVANCING EYECARE