

# MARCO TS-610 WITH FULLY ASSISTED REFRACTION SYSTEM (FARS)

This compact, fully automated system allows patients to conduct the refraction exam themselves. Voice-guided instruction and patient-operated controls allow for self-guided testing that is fast, accurate and reliable...and frees up doctors and technicians to accomplish other tasks simultaneously. **The TRS-610 with FARS will transform your patient workflow, increasing practice efficiencies while providing a more comfortable exam.**

## OPTIMIZED WORKFLOW EFFICIENCIES

Since patients take the test independently, a single staff member can cover multiple devices or multi-task

## TIME-SAVING

Refraction times are greatly reduced, allowing for more patients to be seen per day

## EASY TO USE

Set-up and training is effortless – even an inexperienced technician can facilitate an accurate subjective refraction

## SPACE-SAVING

The compact device mimics a 20' lane in only a 2' x 2' footprint

## IMPROVED PATIENT EXPERIENCE

Patients get to spend less time in the waiting room and testing, leading to more quality time with their doctor



## Fully Assisted Refraction System Specifications

Software	Fully Assisted Refraction System (installation CD)
Joystick controller	
Power supply	5 V DC
Power consumption	100 mA
Dimensions/mass	150 (W) x 100 (D) x 84 (H) mm (including answer lever (upright)) / 0.3 kg 5.9 (W) x 3.9 (D) x 3.3 (H)" / 0.7 lbs.
Voice guidance / Display language	English, French, Spanish, German, Chinese, Japanese
Standard accessories	Velcro tape, Quick reference guide
Compatible device	TS-610
Optional accessories	Control console tray, Tablet stand

## Fully Assisted Refraction System Requirements

Operating system	Windows 10 (32-bit/64-bit) or Windows 11 (64-bit)
Display	Vertical/horizontal resolution: 1,920 x 1,080 pixels or greater
Disk device	CD-ROM drive (only for installation)
Other	Media player installed USB2.0: 1 port or more LAN (100BASE-T or higher): 1 port or more Audio output terminal (if necessary)



### SOLD EXCLUSIVELY BY MARCO

Marco not only offers the best products, but our experts also take the time to understand your vision, providing a comprehensive solution that combines top-notch technology, seamless processes, education, product protection, and expert guidance – that's the Marco difference.



MARCO.COM | 800.874.5274

Marco Healthcare is a division of Advancing Eyecare™

SPEC-MARCO-FARS-V1